Titan BioStar °4

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CE 0124

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Instructions for use

CE 0124

Area of application: Non precious dental disc on pure titanium base, grade 4 for crowns and bridges technique and for metal-ceramic dentures accord. DIN EN ISO 22674, type 4.

Indication: Titan BioStar °4 is a precious metal free, titanium-based alloy for the dental application and used in for single crowns in front and lateral region, bridges in front and lateral region with appropriate cross-section.

Composition (in % by mass)

Ti 99.5

Other constituents:

N,C,H,Fe,O₂

Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.

Technical data (guidelines)

Density g/cm ³	4.5
Elongation at rupture, minimum (%)	18
Tensile strength minimum (MPa):	480
Proof stress 0,2 % (MPa):	483
Modulus of elasticity (Gpa):	120
Melting point (°C):	1.700
WAK (25-500 °C):	9.7 x 10 ⁻⁶ K ⁻¹
Casting temperature (°C)	1.700
Max. burning temperature (°C):	880

Grinding: Grind frames or individual abutments with tungsten carbide burs for titanium. Use also tungsten carbide burs or separating discs for cutting frames off the milling blank. All grinding can be done using cross-cut tungsten carbide burs or diamonds. It is important to cut always in the same direction, so as to avoid metal overlaps that might result in bubbles and porosities when firing the ceramics.

Attention: For top quality and efficiency, please follow all instructions closely. Always use protection gear (goggles, mask, suction...). Titan BioStar °4 dust and chips are highly inflammable! Therefore please follow safety precautions (no smoking when grinding, avoid open flames, etc.). TITAN BioStar °4 is only to be manipulated by trained staff.

Veneering: Use only titanium ceramics. Always follow ceramics manufacturer's instructions. Firing temperature should never exceed 880 °C, since a higher temperature would bring about molecular and structural alterations of the material. After veneering, rubber down all visible metal areas and polish using brushes and polishing pastes.

Veneering: Rubberize parts of framework or crowns or polish to high-gloss with a polishing paste for titanium / titanium based alloys or with rotating tools. In any case follow the instructions of ceramic manufacturer.

Soldering / Welding: To avoid oxide formation during the soldering / welding of titanium only use argon. To protect a corrosion resistance only use additional materials with the same ingredients.

Safety hint: Metal dust is hazardous to health. During the process use personal protective clothing (protecting glasses, dust mask, exhaust). For finishing and sandblasting use a suitable extraction system and/or face mask.

Contraindications, side effects and interactions with other dental alloys: In cases of hypersensitivity (allergy) to the constituents of the alloy, discontinue its use. In individual cases, hypersensitivity reactions (allergies) and electrochemically induced local dysaesthesia have been reported, such as changes in taste and irritation of the oral mucosa. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys. If lasting, electrochemically induced, local dysaesthesia arises from contact with other alloys, the restorations must be replaced with other materials.

Storage conditions: No special storage conditions are required.

Further questions: Contact our technical services (Tel.: +49 (0) 53 21 – 3779-25) or our sales representatives.

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