Zirkon BioStar HT Smile

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CE 0124

gültig für Zirkon BioStar HT Smile / HT Smile ZZ / S-HT Smile / HT Smile Colour / HT Smile Multilayer

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General description: Machinable blank made of yttrium stabilized, presintered zirconium dioxide for milled production of crowns and bridge frameworks on commercial CAD/CAM systems or hand-operated copy-milling machines. Ceramic material according to EN ISO 6872 Typ II, class 4.

Indication: Zirkon BioStar HT Smile is indicated for crowns and bridges (like caps, interlink and all-ceramic restoration with max 3 pontics at anterior and posterior bridges), up to bar constructions, telescopes, Maryland bridges, inlays and onlays.

Contraindications: In the case of bruxism, inadequate hard tooth structure or insufficient space, the use of full ceramic crowns and bridges is not recommended or only advised in some cases. The same applies in the case of inadequate oral hygiene or inadequate coordination of masticatory movements. A full ceramic crown or bridge may also not be used in the event of unsatisfactory preparation results or in the case of production faults.

Preparation information: The minimum wall strength values in the following table are merely reference values. The decision on the manner and undertaking of the preparation is the responsibility of the dental surgeon and dental technician (this includes the validation of the gluing process).

The wall strength should be at least:

Posterior teeth:	occlusal: 0.9 mm	circular 0.5 mm
Anterior teeth:	incisal/occlusal: 0.7 mm	circular 0.5 mm
Abutment crowns of bridge frameworks:	incisal/occlusal 1.0 mm	circular 0.7 mm

The connector surface should be at least:

Anterior bridge with one Pontic:	7 mm ²
Posterior bridge with one Pontic:	9 mm ²

Multilayer blanks are marked with „Incisal/Occlusal“ on the relevant side and have a dent on the bright side.

Milling: Milling can be undertaken with most commercially available CAD/CAM milling machines or hand-operated copy-milling machines. The machine manufacturer's manual should be consulted regarding the selection of the milling cutters and the necessary revolutions/minute. In order to avoid staining or a loss of strength, the blank should if possible not come into contact with liquids (water, oils) or fats (handcream) during milling and related handling of the blank. Wet machining of the dental blanks in their presintered state is in principle possible, but does require that the restoration be dried completely before sintering. Drying can be undertaken in a drying chamber at approx. 80° C or under an infrared lamp. Alternatively, a low-temperature firing program can be selected which will not damage the blank. Notches, sharp edges or slight damage, also on the underside of the framework, may reduce the strength of the completed sintered work. Consequently, it is already important to ensure that surfaces are well rounded and smooth when milling.

Inhalation of dusts may irritate the respiratory tract. Contact with the skin and eyes cause mechanical irritation.

Coloring: If required, the frameworks may be adjusted to the desired tooth color with the standard coloring solutions which are licensed for zirconium oxide ceramics. In order to achieve uniform results, it must be ensured that the framework is completely clean, greaseless and dry before coloring. The instructions of the manufacturers of the coloring solutions should be followed for the further work procedures.

Sintering: Recommended sintering parameters (modifications of these recommended parameters, according to the framework, are possible at your own's authority).

Sintering temperature: 1.450°C ; Soaking time: 2.0 hours ; Heating rate: up to 600 K/h ; Cooling rate: up to 600 K/h

The sintering of the preparations can be undertaken in all standard furnaces which are suitable for zirconium oxide ceramics. The manufacturer's operating instructions must be followed in all cases. The shrinking factor depends on the lot and it is noted on the label, the blank and the inspection certificate. We recommend sintering in a sintering crucible in which the restoration is embedded in suitable granular material. One suitable type of granules is ZrO₂ sintering spheres. When using sintering crucibles and granules, the framework is placed in the granules with the occlusal side downwards. The contact with other frameworks or sinter-bearers should be avoided. If sintering wires or pegs are used, it must be ensured that there is no contact with the cap wall area.

Dental frameworks and crowns made of zirconium oxide should be fired in a furnace which is only used for these products. If additionally restorations made of other ceramic materials are fired or if liquid ceramics are baked in the same furnace, this may impede the sintering behavior or result in localized discoloration. If the firing of different materials in the same furnace cannot be avoided, then an empty firing cycle at maximum temperature is recommended before firing zirconium oxide ceramics.

Facing of the frameworks: Every commercial facing ceramic for frameworks made of zirconium oxide can be used. Please follow the instruction manual of the manufacturer. The CET is about $10.4 \times 10^{-6}/K$.

Finishing after final sintering: Any necessary finishing of the sintered constructions must only be done using water cooling during grinding operations, thus avoiding localized overheating and damage to the material.

Storage information: If the dental blanks are stored in original packaging with regard to temperature and humidity ("cool and dry", namely $10^{\circ}C$ to $35^{\circ}C$ and 30 to 80% relative humidity), then no alteration of the material properties is to be expected, even if stored for longer periods. The blanks must not be stored in direct sunlight.

Disposal of material residues: Remains of the blanks and, of course, also sintered parts made of Zirkon BioStar HT Smile materials are categorized as chemically completely inert and are not subject to any special disposal stipulations.

Further questions: Please contact our technical services (Tel.: +49 (0) 53 21 – 3779-25) or our sales representatives.

Technical Data Sheet

Chemical Formula: $ZrO_2 / Y_2O_3 / Al_2O_3$
Chemical Name: Yttrium stabilized tetragonal zirconium dioxide
Colour: white or coloured, high translucent

		Zirkon BioStar HT Smile (ZZ / S-HT Smile)	Zirkon BioStar HT Smile Colour / HT Smile Multilayer
Chemical characteristics			
ZrO ₂ / HfO ₂ / Y ₂ O ₃	wt%	> 99.0	> 99.0
Y ₂ O ₃	wt%	9.30 ± 0.30	8.50 ± 9.60
HfO ₂	wt%	< 5.0	< 5.0
Al ₂ O ₃	wt%	< 0.1	< 0.1
Fe ₂ O ₃	wt%	< 0.1	< 0.1
Er ₂ O ₃	wt%	none	< 0.1
Density	g/cm ³	> 6.0	> 6.0
4-Point Bending Strength	MPa	600 ± 100	600 ± 100
Thermal expansion	10-6 1/K	10.4	10.4
Open porosity		none	none
Radio activity			
Activity ²³⁸ U	Bq/g	< 0.1	DIN EN ISO 6872
Biocompatibility and chemical solubility			
Cytotoxicity ISO	inoffensive		DIN EN ISO 10993 / DIN EN ISO 9363
Chemical analysis	< 1 µg/cm ²		DIN EN ISO 6872

The data indicated relates to test specimens from which it was obtained and is not unconditionally applicable to other forms of same material. The date must, therefore, be regarded as indicative only.

Date of information: 11/2018

Printing Date: 14.11.2018